





# **Experience the Quality**

## Modern, Clean Design

- Easy to clean due to smooth surfaces, aluminum frame and no slots or rims
- Battery-driven, powerful and silent
- Stationary backrest providing user with a secure feeling
- Patented, self-release suction cups can be easily removed

## Easy to Assemble and Use

- Easy installation without any tools
- Less parts: only 4 main components
- Compact design allows the bath lift to be portable
- Lightweight design

### **Cost-Effective**

Bath lift provides a cost-effective alternative to remodeling bathrooms

## "Experience a Luxurious Bathing Experience with the Invacare® Rio™ Bath Lift."







- 1. Remote Control, Battery Operated
- 2. Textured seating
- 3. Self-Release Suction Cups





# INVACARE® **Rio™ Bath Lift**

It's not just a bath lift...it's a bathing experience.

## **SPECIFICATIONS**

MODEL#	H605
DIMENSIONS	
Seat width without side flaps	14.6"
Seat width with side flaps	28.3"
Seat depth	17.3"
Seat height max.	17.3"
Seat height min.	2.36"
Width of backrest	13.8"-15.7"
Height of backrest	24.2"
Total depth	26.8"
Total height with backrest	40.6"
Base area	26.4" L x 14.6" W
Total weight with battery	21 lb.
Weight, seat	3.5 lb.
Weight, backrest	12.5 lb.
Weight, battery	3.5 lb.
Weight capacity	264 lb.
Color	White

<sup>\*</sup>Each individual should always consult with his/her physician

### Invacare Corporation

www.invacare.com

### USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

© 2011 Invacare
Corporation. All rights
reserved. Trademarks are
identified by the symbols

™ and ®. All trademarks
are owned by or licensed to
Invacare Corporation
unless otherwise noted.
Specifications are subject to
change without notification.

Form No. 11-108 101219



