

**NEW**

Optimum functionality. The Jasmine full body lift has innovative features that provide safer resident handling within a 500 pound weight capacity. Caregivers can rely on safer operation and residents on maximum comfort during everyday transfers.



### Lifting Range

This lift expertly functions in more challenging situations with a rotating 360° 6-point hanger bar and 16.4" minimum height range, essential for assisting residents who may have fallen to the floor.



### Maneuverability

The push handle accommodates any staff member. And with a lighter-weight aluminum and steel construction, control in movement is straightforward and trouble-free.



### Ease For Caregivers

Jasmine lifts feature a battery-powered boom and wide-opening base, both controlled through the hand pendant



### Security

The stable spring-loaded hanger bar produces an audible tone when the sling is engaged. This assures the caregiver that the sling is properly attached.



- Cradle
- XL Cradle



- Scale : JLS5

### Invacare® Jasmine™ Full Body Lift

Height at sling hook-up - MAX	74.6"	Battery	24V DC (Rechargeable)
Height at sling hook-up - MIN	16.4"		2 Shipped with Unit
Base width OPEN	46.7"	Charger input	100-240V AC
Base width CLOSED	25.8"	Charger output/Charging time	29.5 DC Max 6 hours
Base height (Clearance)	4.8"	Audio low battery alarm	Yes
Base length	49"	Motor safety devices	Anti-Entrapment
Caster size (Front/Rear)	4.0"	Limited warranty	3 years on Lift, electronics 2 years on Electronics
Weight capacity	500 lb.	Emergency stop button	Yes
Weight out-of-carton	98 lb.	Emergency up and down	Yes