



Aqua Creek Products
The Total Solution for Pool and Spa Access

The Patriot™ Portable Lift

Item# F-12PPL-HD-A



Inspected And Verified to Meet
Applicable Sections of ADA
Guidelines Section 1009.2 By:

wyle
laboratories

A Registered ISO 9001
Laboratory

The Patriot Portable Pro Pool™ Lift is the first completely portable aquatic lift that is **independently verified to meet ADA guidelines** and capable of lifting a 450lb* (204kg) load.

This lift requires no anchors and can easily be removed from poolside when not needed. It features a compact, low profile frame that does not intrude into the pool when not being used. A unique sliding mechanism allows the lift to be transported easily on its cart and then lowered onto the pool deck for use without using any fastening hardware.

Features:

- 450 Pound* (204 kg) weight capacity
- Independent 3rd party verified to meet ADA standards
- Completely portable unit- no anchor required
- Rechargeable battery with wall mount charger
- Adjustable footrest & lap belt
- Dual flip-up armrests
- Simple push button operation
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty



Options:

- *Upgrade Pack - Item #F-12PPAP*
Includes headrest, chest strap, battery & cover
- *Headrest - Item #F-422PLH*
- *Chest strap - Item #F-423CS*
- *Protective cover - Item #F-12PPC*
- *Pull out leg rest - Item #F-105LAR*
- *Cycle Attachment - Item #F-019CA*



*450 Pound weight capacity for water drafts up to 6"
400 lb. capacity for other applications.

PROUDLY SOLD BY:



Push lift to deck side, set brakes, slide into place and it's ready to use.