Hoyer® Presence

The Hoyer Presence lift has been engineered and designed to handle almost any resident-handling task. Along with the Stature lift, the Presence can lift residents from the floor as well accommodate transfers to higher surfaces. The 500 lbs safe working load and greater spatial area make this lift a perfect fit for larger residents.





Interchangeable Cradle Bars

Part of the restriction of using certain lifts are the slings that you can use on the product. The Presence however, gives you a choice of sling systems.



Powered Base

Operating the Presence is as simple as the push of a button. The base includes two electric actuators which allow the caregiver to open and close the legs without bending or twisting. Tapered leg design allows greater internal area without compromising overall product width.



Hand Pendant

Hand control pads to hold the pendant in place when not in use.



Foot Pad

Ergonomic foot pad to assist with initiating movement.

| Specification | Imperial | Metric |
|------------------------------|----------|---------|
| Safe Working Load | | 227 kgs |
| Maximum Overall Length | 57.5" | 1460 mm |
| Minimum Overall Length | 55.1" | 1400 mm |
| Maximum Overall Height | 81.3" | 2065 mm |
| Minimum Overall Height | 57.5" | 1460 mm |
| Spreader Bar Maximum Height | 75.2" | 1910 mm |
| Spreader Bar Minimum Height | 16.9" | 430 mm |
| Height at Maximum Reach | 51.5" | 1310 mm |
| Reach at Maximum Height | 25.6" | 650 mm |
| Reach at Minimum Height | 15.7" | 400 mm |
| Maximum Reach | 36.6" | 930 mm |
| Turning Radius | | |
| Legs Open - Internal Width | | |
| Legs Closed - External Width | | |
| Front Twin Casters | | |
| Rear Braked Casters | 4.0" | 100 mm |
| Overall Height of Legs | | |
| Ground Clearance | 1.4" | 35 mm |

* Reach = center of spreader bar to the front of the mast

| | 4 Point Position Cradle | (measurement to to | p of location p | oins) |
|--|-------------------------|--------------------|-----------------|-------|
|--|-------------------------|--------------------|-----------------|-------|

| Cradle Max. Height. | | . 61.0" | ' 1550 | mm |
|----------------------|----------|---------|--------|----|
| Cradle Min. Height (| (usable) | . 15.4" | ' 390 | mm |

| Weight | Imperial | Metric |
|--------|----------|--------|
| Total | 88.0 lbs | 40 kgs |

Electric Shock Protection

Charger - Class II

Lift – Internal power source

Degree of Shock Protection

Charger – Type B

Lift – Type B

Intended operating environment: >+5°C <+40°C

Outside this environment functionality and safety may be compromised

Electrical Specifications

Battery - 24vDC rechargeable sealed lead acid type

3.2A hours

Charger rated input - 100-240V AC 24 VDC 50/60Hz Charger rated output - 29.5 VDC, Max. 19 W

