
EZ POOL LIFT^{IM} 10-INCH ANCHOR



THE BOTTOM, AND AN ELECTRO-POLISHED TAMPER RESISTANT LID.

ANCHOR HOUSING: 2-1/2" [7.3 cm OD] SCH-40 304 STAINLESS STEEL PIPE

GROUNDING LUG: BRONZE GROUNDING LUG W/SET SCREW SLEEVE: UHMW (POLYETHYLENE)

LID: ELECTRO-POLISHED 12 GAGE 304 STAINLESS STEEL



9889 GARRY MORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602

*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE *ALL DIMENSIONS ARE IN INCHES

PRODUCT: F-30220		
DWG#: F-30220-1		
SCALE:	1:3	DATE: 12/3/09
DRAWN BY:	DS	REVISION: B