New! MAC (Mobile Aquatic Chair)



Mobile Aquatic Chairs (MAC) are essential for facilities utilizing a ramp, zero-depth entry or movable floor. The MAC is built specifically for use in aquatic environments and with proper maintenance, will look great for years to come.

- · Rear stabilizing wheels provide increased stability when the chair is in the water
- Unibody frame construction for increased durability
- Arm rests flip up easy transfer on and off the chair
- Powder coated stainless steel frame looks beautiful and resists corrosion
- 18" (46 cm) rotomolded plastic seat rigid platform provides transfer support
- 300 pound (136 kg) weight capacity

Model No.	Description	Shipping Wt.	L	W	Н
AC0000	MAC (Mobile Aquatic Chair)	60 lbs.	43"	25"	22"
		(27.2 kg)	(109.5 cm)	(63.5 cm)	(53 cm)









